

DISINTEGRATION TESTER GADT-2B



The GA-DT-2B is designed for Faster and Accurate operation. Operator can test according to USP/EP test protocols. Single or dual baskets can be used. Supplied with complete validation documents (IQ/OQ). Weight of the tablets can be transferred directly from the Balance before and after test for easy calculations of the results.



GALAXY
INTERNATIONAL
+91 40 3580 7244

WWW.GICORP.IN INFO@GICORP.IN
PLOT NO 162, SRI RAMANPURAM STREET NO. 3
UPPAL, HYDERABAD, TELANGANA, INDIA, 500039

SYSTEM SPECIFICATIONS

- **LCD Screen:**
 - Large LCD Screen
- **Keypad:**
 - Multi-Functional Alpha-Numeric Key pad with a smaller number of keys for ease of operation
- **Complies:**
 - USP, IP and other regulatory agencies
- **Features:**
 - Heater protection avoiding frequent failures
 - Provision to operate single or dual baskets
 - Clear visibility of Tablets
 - 50 different methods
- **No of baskets:**
 - 2 baskets of 1000ml each
 - USP type A 2x6 test positions
- **No of strokes:**
 - 30 +/-1 DPM
- **No of strokes:**
 - 55mm +/-2mm
- **Temperature Range:**
 - 5°C to 50°C
 - Accuracy: 0.5°C

ADDITIONAL FEATURES

- External Dot Matrix Printer and Lims connectivity

